

*info@asa-intl.com*



# ASA

*The name of Excellence*


Manufacturer & Exporter All Kinds of

**SURGICAL INSTRUMENTS**



**MHRA**  
Regulating Medicines and Medical Devices

# Wire Extension



**AS-17-171 - 17-177**  
**KIRSCHNER**  
 Wire traction bows,  
 to be used with  
 KIRSCHNER-wires  
 up to 2.5 mm diam.,  
 With 3 hooks,  
 stainless steel

	Inside dimensions width	depth	
17-171	200 mm	155 mm	for the hip
17-172	155 mm	155 mm	for the femur
17-173	135 mm	130 mm	for the knee
17-174	120 mm	120 mm	for the knee
17-175	105 mm	100 mm	for the knee
17-176	95 mm	110 mm	for the calcaneus
17-177	95 mm	70 mm	for the arm and elbow



**AS-17-178**  
 Wrench, double - ended,  
 for wire traction bows



**AS-17-179**  
 Wire tightener, with  
 clamping jaws



17-182	max. 105 mm
17-183	max. 160 mm
17-184	max. 200 mm
17-182	80-95 mm
17-183	90-120 mm
17-184	140-170 mm
diam upto	
17-182	1,6 mm
17-183	2,0 mm
17-184	2,5 mm

**AS-17-182 - 17-184**  
**KIRSCHNER**  
 Wire traction bows, modified,  
 adjustable in width, holds wire  
 securely, with 1 hook



**AS-17-180**  
 Traction hook only,  
 large size




**AS-17-181**  
 Traction hook only,  
 small size



**AS-17-185**  
 Fixation clamp, to secure clamping  
 jaws of 17-179 stainless steel

# Cerclage Instruments




**AS-17-186**  
275 mm

**AS-17-187**  
285 mm

**AS-17-188**  
305 mm



**AS-17-186 - 17-188**  
DEMEL  
Wire guides



**AS-17-189**  
240 mm

**AS-17-190**  
270 mm

**AS-17-191**  
225 mm

**Wire Holding forceps, Wire Tightening Pliers, Flat-nosed Pliers DUROGRIP with Hard Metal Inserts**




**AS-17-192**  
DUROGRIP  
Holding forceps for cerclage wire, with lateral jaws to avoid wire slipping off  
170 mm



**AS-17-193**  
DUROGRIP  
Wire, tightening pliers,  
Ligature pliers  
150 mm



**AS-17-194**  
DUROGRIP  
Wire, tightening pliers,  
Ligature pliers  
150 mm



**AS-17-195**  
with serrated jaws  
140 mm



**AS-17-196**  
one jaw smooth,  
one jaw serrated  
140 mm




**AS-17-197**  
with serrated jaws  
**AS-17-198**  
with longitudinal and  
transversal grooves  
170 mm



**AS-17-199**  
190 mm

# Wire Holding forceps, Flat-nosed Pliers, Pliers for bending and cutting wires




17-200  
130 mm

17-201  
180 mm


AS-17-200 - 17-201  
Extraction pliers  
for boring wires




AS-17-202  
Flat-nosed pliers,  
parallel  
185 mm



AS-17-203  
Vise grip  
200 mm



17-204  
Round Pliers  
180mm



17-204

17-204

17-205

17-206

17-207

17-205 145 mm  
17-206 180 mm  
17-207 180 mm  
Extraction pliers for wires

AS-17-204 - 17-207



AS-17-208  
Bending pliers  
for wires up to  
1.8 mm diam.

# Wire Cutters and wire Cutting Scissors DUROTIP with Carbide Cutting Edges




**AS-17-209**  
ANGLE  
130 mm

for soft wire  
up to diam.  
for hard wire  
up to diam.

0,5 mm  
0,3 mm

**AS-17-210**  
DUROTIP  
ANGLE  
130 mm

0,6 mm  
0,4 mm




**AS-17-211**  
120 mm

**AS-17-212**  
DUROTIP  
120 mm


for soft wire  
up to diam.  
for hard wire  
up to diam.

1,2 mm  
0,8 mm



**AS-17-213**  
HARVEY  
125 mm

0,4 mm  
0,3 mm



**AS-17-214**  
150 mm

for soft wire  
up to diam  
for hard wire  
up to diam

1,5 mm  
1,0 mm



**AS-17-215**  
160 mm

1,5 mm  
1,2 mm



**AS-17-216**  
140 mm  
DUROTIP  
1,5 mm

1,0 mm

# Wire Cutting Pliers, Lateral Cutting Action



**AS-17-217**  
180 mm

for soft wire  
up to diam. 2,0 mm  
For hard wire 1,5 mm  
up to diam  
specially hardened  
cutting edges, front  
and lateral cutting  
action



**AS-17-218**  
REILL  
175 mm

for soft wire  
up to diam. 2,0 mm  
For hard wire 1,6 mm  
up to diam  
specially hardened  
cutting edges, for  
cutting off wires close  
to the bone



**DUROTIP**

**AS-17-219**  
235 mm

for soft wire  
up to diam. 3,0 mm  
For hard wire 2,2 mm  
up to diam

# Wire Cutting Pliers, Lateral and Front Cutting Action




**AS-17-220**  
240 mm

for soft wire up to diam. 2,0 mm  
For hard wire up to diam 1,6 mm



**AS-17-221**  
145 mm

for soft wire up to diam. 1,5 mm  
For hard wire up to diam 1,0 mm



**AS-17-222**  
145 mm

for soft wire up to diam. 2,0 mm  
For hard wire up to diam 1,5 mm



**AS-17-223**  
220 mm

for soft wire up to diam. 3,5 mm  
For hard wire up to diam 2,5 mm




# Wire Cutting Pliers



**AS-17-224**  
320 mm

for soft wire up to diam. 3,0 mm  
For hard wire up to diam. 3,0 mm




Cutter for cutting off KIRSCHNER wires, Tibia head screws, bolts, pins, up to 3 mm diam.



**DUROTIP**

**AS-17-225**  
Cutter for cutting off KIRSCHNER wires, Tibia head screws, bolts, pins, up to 6 mm diam.  
470 mm


# Finger Nail Instruments




**AS-17-226**  
Nail extracting  
forceps  
135 mm




**AS-17-227**  
Nail extracting  
forceps  
145 mm



**AS-17-228**  
Finger nail drill  
with 3 drill points



**AS-17-231**  
Nail splitting  
forceps  
135 mm



**AS-17-232**  
SYSTRUNK  
Scissors, for splitting  
finger nails  
130 mm



**AS-17-233**  
Finger ring saw  
150 mm

# Plaster Shears



AS-17-235  
STILLE-AESCULAP  
260 mm



AS-17-236  
STILLE-AESCULAP  
360 mm

# Plaster Spreader, Plaster Breakers




AS-17-237  
HENING  
270 mm



AS-17-238  
MOD.USA  
230 mm



AS-17-239  
290 mm



with removable jaws



AS-17-240 - 17-241  
WOLFF  
17-240 180 mm  
17-241 240 mm  
plaster breakers



*The name of Excellence*

## **ASA INTERNATIONAL**

*Near Jamia Mosque Usman,  
Kotli Behram Sialkot 51310 Pakistan  
Mob / Whatsapp  
+92 345 6749489 (Swift Response)  
Phone Number  
Tel :+92 52 4267150  
Fax: +92 52 4267151  
Email  
info@asa-intl.net  
Web  
http://asa-intl.net*

## **ASA EXPORT HOUSE**

*Sultan Soap Street,  
Kotli Behram Sialkot 51310 Pakistan  
Mob / Whatsapp  
+92 345 6749489 (Swift Response)  
Phone Number  
Tel :+92 52 4262983  
Fax: +92 52 4267151  
Email  
info@asa-intl.net  
Web  
http://asa-intl.net*



**MHRA**  
Regulating Medicines and Medical Devices

