



# Kerrison Ronguers

Standard, Detachable and Specialty Bone Punches





*ASA International is always committed to introduce excellent quality products in the medical market since 1998. Our efforts to achieve sustainable quality and trustable environment has fetched to the level of great success in Complete Kerrison Range.*

*We Endeavour to serve you more innovated products to the near future.*

*Our Kerrison Ronguers are made of finest stainless steel material which gives 100% performance during work. Please keep a receipt of store from where you bought these Kerrison items and contact them for replacement in case of any defect. Warrantee of this product covers change of Instrument color, hinge problem & breakage can be replaced free of charge.*

*Please note that warrantee does not cover damage or defect caused by misuse, abuse, or failure to use in accordance with the instructions usage provided.*

*Please note that you can only contact our authorize dealers near your area who will contact us to replace your product, if your area is too long in distance then please send defective item by post to authorize dealers, they will guide your replacement parts accordingly to defects and policies.*





KERRISON RONGUERS



Features:

Available in All Sizes

- ☒ 40/90Degree Up biting/down biting
- ☒ Easy Disassembling just in 3 steps without tools
- ☒ Efficient cleaning & Sterilization to increase Instruments life
- ☒ Detailed Marking to indentify Sizes of Kerrisons

Three Easy Steps to Open & Clean Instrument



Step 1:  
Press handle inward

Step 2:  
Press Lock Down

Step 3:  
Leave handle to  
it's maximum opening



Item Number & Tip Size to Identify Instrument size during work



**1. DETACHABLE/NON-DETACHABLE KERRISION RONGUERS (Available in All Sizes)**

**Jaw Position: 40Degree Upwards**

**Detachable/With Lock    Non-Detachable/Without Lock**

Shaft Lengths	Jaw Size Width (mm) Jaw Position 40 Degree Upward	Jaw Opening (mm) Approximate	Standard Footplate	Thin Footplate	Standard Footplate	Thin Footplate	Coating Option (Plasma)
7 Inches	1.0 mm	8.0 mm	ASA-KR-786	ASA-KR-794	ASA-KR-R786	ASA-KR-R794	Black, Blue, Silver
	1.5 mm	9.0 mm	ASA-KR-787	ASA-KR-795	ASA-KR-R787	ASA-KR-R795	Black, Blue, Silver
	2.0 mm	9.0 mm	ASA-KR-788	ASA-KR-796	ASA-KR-R788	ASA-KR-R796	Black, Blue, Silver
	2.5 mm	10.0 mm	ASA-KR-789	ASA-KR-797	ASA-KR-R789	ASA-KR-R797	Black, Blue, Silver
	3.0 mm	10.0 mm	ASA-KR-790	ASA-KR-798	ASA-KR-R790	ASA-KR-R798	Black, Blue, Silver
	4.0 mm	12.0 mm	ASA-KR-791	ASA-KR-799	ASA-KR-R791	ASA-KR-R799	Black, Blue, Silver
	5.0 mm	12.0 mm	ASA-KR-792	ASA-KR-800	ASA-KR-R792	ASA-KR-R800	Black, Blue, Silver
	6.0 mm	12.0 mm	ASA-KR-793	ASA-KR-801	ASA-KR-R793	ASA-KR-R801	Black, Blue, Silver

**Jaw Position: 40Degree Upwards**

**Detachable/With Lock    Non-Detachable/Without Lock**

Shaft Lengths	Jaw Size Width (mm) Jaw Position 40 Degree Upward	Jaw Opening (mm) Approximate	Standard Footplate	Thin Footplate	Standard Footplate	Thin Footplate	Coating Option (Plasma)
8 Inches	1.0 mm	8.0 mm	ASA-KR-802	ASA-KR-810	ASA-KR-R802	ASA-KR-R810	Black, Blue, Silver
	1.5 mm	9.0 mm	ASA-KR-803	ASA-KR-811	ASA-KR-R803	ASA-KR-R811	Black, Blue, Silver
	2.0 mm	9.0 mm	ASA-KR-804	ASA-KR-812	ASA-KR-R804	ASA-KR-R812	Black, Blue, Silver
	2.5 mm	10.0 mm	ASA-KR-805	ASA-KR-813	ASA-KR-R805	ASA-KR-R813	Black, Blue, Silver
	3.0 mm	10.0 mm	ASA-KR-806	ASA-KR-814	ASA-KR-R806	ASA-KR-R814	Black, Blue, Silver
	4.0 mm	12.0 mm	ASA-KR-807	ASA-KR-815	ASA-KR-R807	ASA-KR-R815	Black, Blue, Silver
	5.0 mm	12.0 mm	ASA-KR-808	ASA-KR-816	ASA-KR-R808	ASA-KR-R816	Black, Blue, Silver
	6.0 mm	12.0 mm	ASA-KR-809	ASA-KR-817	ASA-KR-R809	ASA-KR-R817	Black, Blue, Silver



**Jaw Position: 40Degree Upwards**

**Detachable/With Lock   Non-Detachable/Without Lock**

Shaft Lengths	Jaw Size Width (mm) Jaw Position 40 Degree Upward	Jaw Opening (mm) Approximate	Standard Footplate	Thin Footplate	Standard Footplate	Thin Footplate	Coating Option (Plasma)
9 Inches	1.0 mm	8.0 mm	ASA-KR-818	ASA-KR-823	ASA-KR-R818	ASA-KR-R823	Black, Blue, Silver
	2.0 mm	9.0 mm	ASA-KR-819	ASA-KR-824	ASA-KR-R819	ASA-KR-R824	Black, Blue, Silver
	3.0 mm	10.0 mm	ASA-KR-820	ASA-KR-825	ASA-KR-R820	ASA-KR-R825	Black, Blue, Silver
	4.0 mm	12.0 mm	ASA-KR-821	ASA-KR-826	ASA-KR-R821	ASA-KR-R826	Black, Blue, Silver
	5.0 mm	12.0 mm	ASA-KR-822	ASA-KR-827	ASA-KR-R822	ASA-KR-R827	Black, Blue, Silver

**Jaw Position: 40Degree Upwards**

**Detachable/With Lock   Non-Detachable/Without Lock**

Shaft Lengths	Jaw Size Width (mm) Jaw Position 40 Degree Upward	Jaw Opening (mm) Approximate	Standard Footplate	Thin Footplate	Standard Footplate	Thin Footplate	Coating Option (Plasma)
11 Inches	1.0 mm	8.0 mm	ASA-KR-828	ASA-KR-833	ASA-KR-R828	ASA-KR-R833	Black, Blue, Silver
	2.0 mm	9.0 mm	ASA-KR-829	ASA-KR-834	ASA-KR-R829	ASA-KR-R834	Black, Blue, Silver
	3.0 mm	10.0 mm	ASA-KR-830	ASA-KR-835	ASA-KR-R830	ASA-KR-R835	Black, Blue, Silver
	4.0 mm	12.0 mm	ASA-KR-831	ASA-KR-836	ASA-KR-R831	ASA-KR-R836	Black, Blue, Silver
	5.0 mm	12.0 mm	ASA-KR-832	ASA-KR-837	ASA-KR-R832	ASA-KR-R837	Black, Blue, Silver



Jaw Position: 40Degree Downwards

Detachable/With Lock    Non-Detachable/Without Lock

Shaft Lengths	Jaw Size Width (mm) Jaw Position 40 Degree Upward	Jaw Opening (mm) Approximate	Standard Footplate	Thin Footplate	Standard Footplate	Thin Footplate	Coating Option (Plasma)
7 Inches	1.0 mm	8.0 mm	ASA-KR-838	ASA-KR-843	ASA-KR-R838	ASA-KR-R843	Black, Blue, Silver
	2.0 mm	9.0 mm	ASA-KR-839	ASA-KR-844	ASA-KR-R839	ASA-KR-R844	Black, Blue, Silver
	3.0 mm	10.0 mm	ASA-KR-840	ASA-KR-845	ASA-KR-R840	ASA-KR-R845	Black, Blue, Silver
	4.0 mm	12.0 mm	ASA-KR-841	ASA-KR-846	ASA-KR-R841	ASA-KR-R846	Black, Blue, Silver
	5.0 mm	12.0 mm	ASA-KR-842	ASA-KR-847	ASA-KR-R842	ASA-KR-R847	Black, Blue, Silver

Jaw Position: 40Degree Downwards

Non-Detachable/Without Lock

Shaft Lengths	Jaw Size Width (mm) Jaw Position 40 Degree Upward	Jaw Opening (mm) Approximate	Standard Footplate	Thin Footplate	Coating Option (Plasma)
9 Inches	2.0 mm	8.0 mm	ASA-KR-848		Black, Blue, Silver
	3.0 mm	10.0 mm	ASA-KR-849		Black, Blue, Silver
	4.0 mm	12.0 mm	ASA-KR-850	NIL	Black, Blue, Silver
	5.0 mm	12.0 mm	ASA-KR-851		Black, Blue, Silver

